

## COMPATIBILITY BETWEEN PEN-NEEDLE AND PEN INJECTOR DEVICES

To whom it may concern, this statement certifies that Insupen® Original and Advanced range of pen needles, which are manufactured in Italy in accordance with EN ISO 11608-2:2022, are compliant to the standard above mentioned for compatibility between pen needle and needle injection system (paragraph 9.4).

This statement means that Insupen® Original and Advanced can be used in combination with the following injection devices:

<p>AstraZeneca/Amylin Pharm</p> <ul style="list-style-type: none"> <li>• Byetta Pen 5 mcg</li> <li>• Byetta Pen 10 mcg</li> <li>• Symlin Pen 60</li> <li>• Symlin Pen 120</li> </ul> <p>Berlin Chemie:</p> <ul style="list-style-type: none"> <li>• BerliPen Areo 2</li> <li>• BerlinPen Areo 3</li> </ul> <p>BD:</p> <ul style="list-style-type: none"> <li>• Wangbangpen</li> </ul> <p>Bioton:</p> <ul style="list-style-type: none"> <li>• GensuPen</li> <li>• GensuPen Improve</li> <li>• Genupen 2</li> </ul> <p>Brighter</p> <ul style="list-style-type: none"> <li>• Actiste</li> </ul> <p>Eli Lilly:</p> <ul style="list-style-type: none"> <li>• HumaPen Luxura</li> <li>• HumaPen Luxura HD</li> <li>• HumaPen Memoir</li> <li>• Humalog KwikPen</li> <li>• Humalog Mix 50/50 KwikPen</li> <li>• Humalog Mix 75/25 KwikPen</li> <li>• Humalog U200 KwikPen</li> <li>• Basaglar KwikPen</li> <li>• HumaPen Savvio</li> <li>• Forsteo</li> <li>• Humulin N KwikPen</li> <li>• Humulin R KwikPen</li> <li>• Humulin 70/30 KwikPen</li> </ul> <p>Owen Mumford:</p> <ul style="list-style-type: none"> <li>• AutoPen Classic 1-21</li> <li>• Autopen Classic 2-42</li> <li>• AutoPen 24 1-21</li> <li>• AutoPen 24 2-42</li> </ul>	<p>Ypsomed:</p> <ul style="list-style-type: none"> <li>• Ypsopen</li> </ul> <p>Mylan:</p> <ul style="list-style-type: none"> <li>• Semglee</li> </ul> <p>Novo Nordisk:</p> <ul style="list-style-type: none"> <li>• NovoPen 3</li> <li>• NovoPen 4</li> <li>• NovoPen 5</li> <li>• NovoPen Echo</li> <li>• NovoPen Junior</li> <li>• Novo Rapid Flexpen</li> <li>• Levemir Flexpen</li> <li>• Novolog Flexpen</li> <li>• Novolog Mix 70/30 Flexpen</li> <li>• Innolet</li> <li>• Victoza</li> <li>• Flex Touch</li> <li>• FlexTouch Tresiba Pen U100</li> <li>• Flex Touch Tresiba Pen U200</li> </ul> <p>Sanofi:</p> <ul style="list-style-type: none"> <li>• AllStar</li> <li>• AllStar Pro</li> <li>• KlikStar</li> <li>• JuniorStar</li> <li>• Lyxumia Adlyxin Pen</li> <li>• Suliqua 10-40 SoloStar</li> <li>• Suliqua 30-60 SoloStar</li> <li>• Soliqua 100/33</li> <li>• Lantus SoloStar</li> <li>• Apidra SoloStar</li> <li>• Insuman Basal SoloStar</li> <li>• Insuman Rapid SoloStar</li> <li>• Insuman Comb25 SoloStar</li> <li>• Toujeo SoloStar</li> <li>• Toujeo Max SoloStar</li> <li>• Toujeo DoubleStar</li> <li>• Insulin lispro SoloStar</li> <li>• Admelog SoloStar</li> <li>• Insulin Aspart</li> <li>• TactiPen</li> </ul> <p>Tonghua Dongbao:</p> <ul style="list-style-type: none"> <li>• Gansulin Pen</li> </ul>
---	---

All data are held on file by Pikdare S.p.A. to support this verification.

Date: 06.03.2023

**HTL-STREFA, Inc.**

113 Towne Lake Pkwy Ste 120, Woodstock, GA 30188, USA Tel.: +1 770 528 0410

[www.mtdglobal.com](http://www.mtdglobal.com) [www.mydroplet.com](http://www.mydroplet.com) [www.dropsafe.info](http://www.dropsafe.info)

[info@htl-strefa.com](mailto:info@htl-strefa.com)

**Medical Technology and Devices S.p.A.**

Registered office: Via Saldarini Catelli, 10 – 22070 Casnate con Bernate (CO) – Italy

Tel: 0039 031 7297111 – [www.mtdglobal.com](http://www.mtdglobal.com)