

Droplet pen needles fit the major brands of pen injectors in the US and Canada*

droplet®

Insulin

Novo Nordisk



FlexPen



Flex Touch



NovoPen®



NovoPen® Junior



Tresiba® FlexTouch®



Levemir® FlexPen®



Innolet



NovoLog® FlexPen®



NovoRapid® FlexPen®



Xultophy

Eli Lilly



HumaPen®



Humulin Kwik Pen



Forsteo



Humalog® KwikPen®



Humalog®
U-200 KwikPen®



Humalog®
Mix 50/50™ Kwik Pen®



Humalog®
Mix 75/25™ Kwik Pen®

Sanofi-Aventis



AllStar®



SoloSTAR®



Lyxumia Pen



ClikSTAR®



OptiPen®

B. Braun



Omnican® Pen 31



Omnican® Pen 32

Haselmeier



BerliPen® JUNIOR



BerliPen®



JuniorSTAR



TactiPen®



i-pen

Owen Mumford



BerliPen® 301 & 302



Autopen®

Ypsomed



Ypsopen®

GLP-1

Novo Nordisk



Ozempic® 1mg



Saxenda® 6mg/ml



Xultophy



Wegovy® 2,4mg FlexTouch®



Victoza Pen

Eli Lilly



Mounjaro® 15mg/dose
KwikPen®



Zepbound Kwikpen

* For a complete compatibility list, scan here.



Droplet pen needles are designed to create a smoother, more comfortable injection experience.

droplet[®]



Micron
34G x 3.5mm (9/64")
100 pc/box



32G x 8mm (5/16")
100 pc/box



32G x 4mm (5/32")
100 pc/box



31G x 5mm (3/16")
100 pc/box



32G x 5mm (3/16")
100 pc/box



31G x 6mm (1/4")
100 pc/box



32G x 6mm (1/4")
100 pc/box



31G x 8mm (5/16")
100 pc/box



29G x 12mm (1/2")
100 pc/box

**Feels like nothing.
Delivers like nothing else.**